




Pro-gard Products

DTK™

Procedure Card for TP1001


General Drug Screening**Read Before Use**

As a precaution, gloves and eye protection should be worn when handling drug substances and using test units. Do not puncture any part of unit. Do not remove probe once a capsule is broken. In the unlikely event of a spill, rinse affected area thoroughly with water. Keep out of reach of children. All results are presumptive.


Test Preparation**1** Remove the probe from test unit.**2** Sample suspected substance or surface.

Touch probe to the suspected substance.

Or swipe either side of the probe across the suspect surface.

**3** Secure probe into test unit. Press firmly

Do not continue unless probe is secured properly in place.

Chamber 1 Identification**4** From both sides of the kit, squeeze on the circle of **chamber 1** until the capsule breaks.**5** Gently shake test unit for 2-3 seconds.**Chamber 1 Results**


If a brown powder is present, the appearance of a purple-violet color strongly indicates the presence of Brown Heroin. If a brown powder is not present, the appearance of a purple-violet color indicates either the presence of basic alkaloid substance (possible drug - see below).



The appearance of a green or blue-green color indicates the presence of **Crack**.



If any of the above colors develop, proceed to Disposal. If none of the above colors develop, proceed to Step 6.

Chamber 2 Identification**6** From both sides of the kit, squeeze on the circle of **chamber 2** until the capsule breaks.**7** Gently shake test unit for 2-3 seconds.**Chamber 2 Results**

The appearance of a purple-violet color indicates the presence of one of these suspected drugs:

Heroin, Amphetamine, Methamphetamine, Ecstasy/MDMA, Buphedrone, Cathinone, DMT, Fentanyl, Ketamine, Mephedrone, Methadone, MethCathinone, Methylone, mCPP, MDPBP, MDPV, MPA, PCP, PMA, PMMA, 2-aminoindane (2-AI), other synthetic/designer drugs.



The appearance of a green or blue-green color indicates the presence of **Cocaine**.



If none of the above colors develop, none of the suspected drugs listed above are present.

**Disposal**

The test kit is safe to dispose of in the trash.

04-21